




Technical developments : RFQ/Laser





TITAN RFQ – A digital linear Paul trap

Quick facts and figures

- 700 mm long, $r_0 = 10$ mm
- $C \approx 1500$ pF
- Stack of optically triggered MOSFETs produces RF
- 200 kHz to 1200 kHz frequency range
- Amplitude up to $800 V_{PP}$



Gas-filled (He,H) linear Paul trap



TITAN RFQ – Forward extraction

TITAN RFQ primarily needed to:

- Decelerate ISAC's beam from < 40 keV to 2 keV
- Cool the incoming beam (reduce the phase space volume)
- Bunch the incoming DC beam and send pulses to MPET (EBIT)

Pulsed drift tube

- Defines beam energy
- Switches ions to GND potential

Incoming
beam energy


20 keV



RFQ
20 kV



PDT
18 kV





Ion
elevator

↑
V

GND



Outgoing beam
energy 2 keV



Status (T. Brunner SMI '10)

- Fully operational at 20 kV (^8He beam time with 3 ions/minute at MPET MCP)
- Commissioned for 40 kV
- Frequency range from 250 kHz to 1200 kHz
- DC transmission of up to 80 % for Cs
- Broad mass range demonstrated for ion masses from 6 to 133
- Cooling with He and H possible
- Several online beam times with radioactive He, Li, K, Rb, Ca, In, Cs

Nim B. in preparation (T. Brunner et al.)

This is as far as forward extraction goes....

Status (T. Brunner SMI '10)

- Fully operational at 20 kV (^8He beam time with 3 ions/minute at MPET MCP)
- Commissioned for 40 kV
- Frequency range from 250 kHz to 1200 kHz
- DC transmission of up to 80 % for Cs
- Broad mass range demonstrated for ion masses from 6 to 133
- Cooling with He and H possible
- Several online beam times with radioactive He, Li, K, Rb, Ca, In, Cs
Nim B. in preparation (T. Brunner et al.)

This is as far as forward extraction goes....




...this RFQ is also able to reverse-extract beams (demonstrated circa '06)

Which is very handy for laser spectroscopy :)

Reverse-bunched-beam laser spectroscopy

1) Cool, trap and reverse extract bunch



2) Only count when the bunch arrives



Background reduced by



$$\frac{t_{acc}}{t_{bunch}} \approx \frac{10ms}{1.6\mu s} \approx 6.25 \times 10^3$$

Fluorescence detection




~ 10⁵ ions (50 Hz cycle)


First on-line data on ^{87}Rb ($I=0,4$), $\sim 1\text{pA}$



First on-line data on ^{78,78m}Rb (I=0,4), ~ 1pA



Laser : 0.35 mW and Disp. ~ 6 MHz/V)
Lorentz = 7.6(1.6) /Gauss = 28.5(2.1)



Hyperfine structure: $F = I + J$

$$W_F = W_J + A \frac{C}{2} + B \frac{(3C/4)(C+1) - I(I+1)J(J+1)}{2I(2I+1)J(2J+1)}$$

$$A = \frac{\mu_I B_{el}}{I J} \quad B = e Q_s \left\langle \frac{\partial^2 V_e}{\partial z^2} \right\rangle$$

$$C = F(F+1) - I(I+1) - J(J+1)$$


Results (in MHz)

	IS (⁷⁸ - ^{78m})	A (S _{1/2})	A (P _{3/2})	B (P _{3/2})
This work	67.1 (2.1)	1184.4 (0.6)	29.5 (0.2)	77.1 (2.2)
Other work ¹	74.6 (2.2)	1186.8 (1.0)	29.4 (0.5)	76.5 (3.7)

[1] C. Thibault et al. PRC 23 6 (1981)

This is not enough (never is)...so we need to make improvements


Ars longa,
vita brevis,
occasio praeceps,
experimentum periculosum.
iudicium difficile.




(Hippocrates)

A closer look at the polarizer line


P. Levy improved the design > lifetime



A closer look at the polarizer line





A closer look at the polarizer line





Seeing the laser mode at the interaction region

Coherent CCD camera + beam analyser resurrected:

- Test made with 670.9 nm (camera is ok for 780 nm)
- > 90 % of the laser mode within in 1mm




TITAN switch yard - 100 Hz rep. rate





Switching: +3.5kV to -3.5kV

Glassman 7.5mA Supply, 25kOhm Resistor, Short cable




- > 1ms settling time achieved*
- 4 x Glassman 7.5 mA supp required
- Integration with EPICS
- 20/80 duty cycle




* Supplies 1.5 m away from switch box

Another look at the polarizer line



- Doppler shift for ^{78}Rb @ 28 KeV $\sim 0.7 \text{ nm}$
-> Background reduced only by a factor of 2...
- Might also be throwing away signal... (ang. dependence of interference filter)

Another look at the polarizer line



- **Less background and no need of interference filter**
- **...but solution requires fibre coupling!!**

My wish list

- **Measure the beam emittance (longitudinal and transverse) at the polarizer line before and after RFQ**
- **Voltage stability of RFQ (measure hfs over long t)**
- **Space-charge limit**
- **Optimal settings for operation (cooling time, RF, etc)**
- **Investigate double trap issue...**


For that we need OLIS

OLIS beam time essential for the laser programme :

- Find tuning that works best to maximise transport and laser /ion overlap
- **Efficiency measurements for Rb experiment (laser co and counter propagating)**

A closer look at the 78Rb spectra

Ground state



Strongest peak isomer

